Laminar Flow Cabinet Vd-650



Ningbo, China

Port:Material:Stainless SteelType:Slit TypeFunction:ExhaustFeature:Acid & Alkali ResistantHood Type:StandardColor:Metal**Gontact NowHease Contact NotBasic Info.**Metal

Model NO.
VD-650
Customized
Non-Customized
Condition
New
Inlet Pressure
>0.5PA
>0.5PA

Trademark
SUGOLD
Transport Package
Wooden
Specification
640*530*845mm
Origin
Zhejiang, China
HS Code
8414809053

Product Description

laminar flow cabinet VD-650

Vertical air flow clean bench

It is widely used in medicine, pharmacy, biological & chemical experiment, electron, precise apparatus, instrument and so on.


Contact switch adjusts voltage to keep wind speed in working area in ideal state all the time.

Technical Parameters				
Parameter\Model	VD-650	VD-850	HD-650	HD-850
Clean grade	Grade 100(209E U.S. Federal)			
Number of bacteria	<0.5per utensil.hour(90mm utensil)			
Noise	<62dB			
Powersupply	AC,220V/50Hz			
Equiblast Speed	0.25-0.45m/s			
Dimension of working area(W1×D1×H1)(mm)	485×475×455	685×475×455	485×385×455	685×385×455
Overall dimension(W1×D1×H1)(mm)	640×530×845	840×530×845	640×530×845	840×530×845
Specification and number of light from firefly/ ultraviolet light	8W×1/10W×1	8W×1/10W×1	8W×1/10W×1	8W×1/10W×1
Specification and number of high effective filter	480×460×38	680×460×38×1	480×480×38×1	680×480×38×1

Note: Special specification of non-standard type clean bench is also available on demand.







Horizontal air supply