WiFi Intraoral Camera Endoscope Led Light USB Micro Cam Dental Teeth Photo Shoot Video Home Use Oral Cam Dentist IP67 Waterproof



Wireless Oral Dental Wifi Intraoral Camera Endoscope Borescope Led Light USB Camera Teeth Photo Shoot Dentist Intra Oral Camera

Features:

- 1) Simple, Light, Streamlined appearance.
- 2) Camera can be rotated at 360 degree.
- 3) Split design, easy to clean and disinfect.
- 4) Waterproof, Dustproof, High brightness LED lights.
- 5) It can take pictures and videos, and save in the Mobile phone or Computer.

6) Free Apps for Android/iOS/Tablet/Windows system; support 4 users at the same time.

Parameter:

WIFI distance: 10 meters Battery Life: About 2 hours Taking pictures and Video: Support Power: 3.7V Resolution: 1280*720 , 640x480 Pixel: 1.3MP Auxiliary light source: 6pcs LED bright lights Lens: F / NO 2.8 Visual angle: 100 degree Magnification: 150x Waterproof Level: IP67 Charging time: About 1 hour Sign of fully charged: 4pcs Blue LED lights All ON Supported OS: Android 4.4 and above, iOS 9.0 and above, Windows Weight: 170g

Package List:

- 1 x oral endoscope (Camera + Handle)
- 1 x charging cable
- 1 x User manual

Note: Comes With or Without Retail Package (Optional)

Tips:

- 1. When only one power indicator light is bright, please give it charge.
- 2. Please use a 5V/1A charger to charge.
- 3. While charging, the blue Indicator light flashes.
- 4. When all four blue power indicators are lit, indicates that the full charge is completed.





